DRUG ELUTING BALLOON CATHETER

AcoArt Orchid
Drug Eluting Balloon 0.035"

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DEVICE NAME


The drug eluting balloon catheter (DEB): AcoArt Orchid (DEB 0.035"), AcoArt Tulip (DEB 0.018") and AcoArt Litos (DEB 0.014") are percutaneous transluminal angioplasty (PTA) catheters with different guide-wire compatibilities and are manufactured by Acotec Scientific Co., Ltd.

DESCRIPTION

The DEB catheters are over the wire (OTW) peripheral balloon catheter, specifically designed for percutaneous transluminal angioplasty in atherosclerotic obstructed vessels. The catheter has a dual-lumen shaft, onto the distal tip of which a balloon is welded. The double-lumen shaft is branched at the proximal end so that one tube forms the entrance to the central lumen for the guide wire, while the other tube is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution.

DEB catheters are paclitaxel coated balloon dilatation catheter, the dose of paclitaxel is 3 µg per square millimeter.

The DEB catheters are available in different balloon sizes. Nominal balloon diameter/ lengths and guide wire compatibility are printed on the hub.

INDICATIONS

The drug eluting balloon (DEB) catheters are indicated for percutaneous transluminal angioplasty (PTA) in patients with obstructive disease of peripheral arteries.

CONTRAINDICATIONS

- The DEB catheters are contraindicated for the use in coronary arteries, supra-aortic vessels and cerebrovascular.
- Inability to cross lesion with a guide wire.
- The DEB catheters should not be used in pregnant or breast feeding women or in patients with known hypersensitivity to paclitaxel.

WARNINGS

- This device is designed and intended for single use only. **DO NOT RSTERILIZE AND/OR REUSE.** Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. ACOTEC will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse. The coating excludes any reuse of the device.
- **The maximal treated lesion length should not be more than 75 cm to avoid potential systemic effects of paclitaxel (e.g. myelosuppression). If the sum of lesions is longer than 75 cm a second intervention has to be scheduled at a later point in time.**
- Inspect the device, prior to procedure, to verify functionality and lack of damaged parts. Do not use the device if the outer or the inner package is damaged or opened.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just distal to the stenosis.
- When the catheter is in the body, it should be manipulated while under sufficient and/or high quality fluoroscopy. Prior to withdrawing the balloon catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Never use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium.
- Do not expose the device to organic solvents, e.g. alcohol
- Do not manipulate the balloon catheter in inflated state. The position of the balloon catheter may only be changed with the guide wire in place.
- If resistance occurs during manipulation, the cause must first be ascertained by fluoroscopy, road mapping or DSA before the balloon catheter is moved backwards or forwards.
- The guide wire may under no circumstances be moved during inflation of the balloon.
- The balloon must be completely deflated before retrieving the catheter from the vascular system.

- Do not exceed the Rated Burst Pressure (RBP). The RBP is based on the results of in-vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the rated burst pressure may cause the balloon to rupture.
- Only the physician trained and after training of the balloon dilatation catheter could use the device.
- To reduce the risk of vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just distal to the stenosis.
- Use only a mixture of contrast medium and saline solution to fill the balloon (1:1). Never use air or any gaseous medium to inflate the balloon of the catheter.
- Use by the labeled date.
- Do not use with Lipiodol™ and Ethiodol™ contrast media, or similar relevant ingredients contrast agents.
- Never apply positive pressure to the balloon during preparation.

PRECAUTIONS

- Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for PTA before insertion of the dilatation catheter.
- To minimize the possible introduction of air into the system, it is imperative that before proceeding careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Precautions to prevent or reduce clotting should be taken when any catheter is used. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.
- The balloon catheter should be used with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- Allergic reactions to contrast medium should be identified before treatment.
- Catheter applications vary and the technique must be selected on the basis of the patient's condition and the experience of the interventionalist.
- Never advance the balloon catheter without the guide wire extending from the tip.
- Never attempt to move the guide wire when the balloon is inflated.
- Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- Store at controlled room temperature, in a dry place. Keep away from sunlight.
- Don't exceed dilatation pressures to "rated burst pressure". In any case, always use manometer controlled inflation.

CAUTION: Larger sizes of the balloon catheter may exhibit slower deflation times particularly on long catheter shafts.

POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

Complications associated with the use of the DEB balloon catheter are similar to the ones associated with standard PTA procedures. Possible complications may include, but are not limited to:

Puncture related:

- Local hematoma
- Local hemorrhage
- Local or distal thromboembolic episodes
- Thrombosis
- Arterio-venous fistula
- Pseudoaneurysm
- Local infections

Dilatation related

- Dissection in the dilated artery wall
- Perforation of the artery wall
- Prolonged spasms
- Acute re-occlusion necessitating surgical intervention
- Restenosis of the dilated artery
- Total occlusion of the peripheral artery

Angiography related

- Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection

- Systemic embolization
- Endocarditis
- Short-term hemodynamic deterioration
- Death
- Drug reactions
- Allergic reaction to contrast medium
- Pyrogenic reaction

PHARMACOLOGICAL INTERACTION

Due to the low dose and local administration, the pharmacological interactions are not to be expected

INSTRUCTIONS FOR USE

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the dilatation catheter, to verify proper function. Verify that the catheter and sterile packaging have not been damaged in shipment and that the catheter size is suitable for the specific procedure for which it is intended.

Handle the device with extreme caution in order to avoid any damage to the folded balloon. Moistening of the device with saline, insertion through the sheath and contact with blood will not reduce the efficacy of the drug coating as long as the balloon remains uninflated.

1. PREPARATION EXPANDING DEVICE

- Prepare the inflation device according to the manufacturer's instructions.
- Carefully remove the balloon protective sheath. Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment. In case of resistance, twist the protection with one hand holding the shaft with the other.

2. CATHETER SELECTION

- The nominal balloon size must be chosen equal to the inner diameter of the artery distal to the lesion. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter (non drug coated) catheter to pre-dilate the lesion to facilitate the passage of a more appropriate-sized dilatation catheter.
- In cases of total occlusion or subocclusive lesions a predilatation with a (non drug coated) balloon is recommended.

3. DILATATION CATHETER PREPARATION

- The catheter is packaged in a protective tube. Carefully remove the catheter from the protective tube.
- The folded balloon catheter contains tiny air bubbles that should be purged prior use. To do this, keep the balloon catheter in a downward vertical position. Keep the protective sheath in place during the purging procedure. Connect a Luer-lock syringe partially filled with saline-contrast mixture solution to the inflation port of the catheter hub. Apply negative pressure for 15 seconds until air is completely evacuated and release the plunger. Repeat this operation a couple of times until migration of air bubbles towards the syringe stops. **Never use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium. Never apply positive pressure to the balloon during preparation.**
- Remove the protective tube from the balloon.
- Flush the wire lumen properly through the connector luer connector.

4. INFLATION DEVICE CONNECTION TO CATHETER

- To remove any air lodged in the distal Luer fitting of the inflation device, purge approximately 1 ml (cc) of contrast medium.
- With the stopcock in the closed position, disconnect the syringe used in preparation applying a slight positive pressure. A meniscus of contrast medium will appear in the balloon port when the syringe is removed. Verify that a meniscus of contrast medium is evident in both the dilatation catheter balloon port (hub) and the inflation device connection. Securely couple the inflation device to the balloon port of the balloon dilatation catheter.

5. USE OF THE DRUG ELUTING BALLOON CATHETER

- Insert a guidewire through the hemostatic valve following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guiding catheter/ introducer sheath. When complete, withdraw the guidewire/ introducer, if used.
- Attach a torque device to the wire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis.
- Thoroughly aspirate and flush the guiding catheter in preparation (when guiding catheter is used) for introduction of the dilatation catheter. Back load the distal tip of the dilatation catheter onto the guidewire.

NOTE: To avoid kinking, advance the dilatation catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.

- Advance the catheter through the hemostatic valve slowly, while the balloon is fully deflated. It should be observed that the hemostatic valve is only closed as much as to prevent blood return yet permitting easy movements of the dilatation catheter. If resistance is encountered, do not advance the catheter through the adapter.
 - Under fluoroscopy, use the balloon radiopaque markers to position the balloon within the lesion to be dilated and inflate the balloon to the appropriate pressure (refer to balloon compliance table).
 - The majority of the drug is released within the first 30 seconds of balloon inflation. The duration of the inflation should therefore be between 30 seconds and 1 minute for optimal drug release. In order to optimize lesion dilatation, longer inflation times are possible at the discretion of the operator.
 - Completely deflate the balloon catheter. Withdraw the deflated dilatation catheter and guidewire from the guiding catheter/introducer sheath, through the hemostatic valve. Tighten the knurled knob on the hemostatic valve.
 - If necessary the balloon catheter can be exchanged over the guidewire which remains in the vessel, for different balloon types or sizes.
- WARNING:** The same vessel segment must not be dilated with more than one drug eluting balloon. Should a postdilatation be necessary after the utilization of a drug eluting balloon, so should this be done with a standard (non drug coated) dilatation balloon. Drug Eluting Stents must not be implanted into the vessel segment that has been treated with a drug eluting balloon.
- If more than one balloon is required to treat a single lesion, balloons must overlap by at least 1 cm. A second new balloon catheter has to be used because the drug is almost completely released during the first expansion. Additional overlapping should be avoided.

6. PRE- AND POST- PROCEDURE ANTIPLATELET REGIMEN

- Dual antiplatelet therapy (ASS and clopidogrel or ticlopidine) should be administered pre-procedure and for a minimum of 4 weeks after the intervention. Prolonged antiplatelet therapy can be given at the discretion of the physician and should be considered after placement of stents.

PRODUCT SUPPLY

The DEB catheters are supplied sterile and intended for single use only. The DEB catheters are sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use product prior to labeled Expiration Date.

CAUTION: Do not use if the inner package is open or damaged.

SHIPMENT CONDITION

Product should be shipped, avoiding heavy weight, direct sunlight, snow and rain.

STORAGE CONDITION

The packed catheter shall be stored in a dry, cool and clean place free of corrosive gases.

WARRANTY/LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in Acotec's instructions for use are expressly to be considered as an integral part of this provision. Acotec guarantees the product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Acotec disclaims any warranty of merchantability or fitness for a particular purpose of the product. Acotec is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on Acotec's part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Acotec, however originating, from the product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Acotec's legitimate interest in limiting its liability or warranty. No person has any authority to bind Acotec to any warranty or liability regarding the product.

Additional information regarding paclitaxel medical devices Warning

A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.

Summary of the meta-analysis

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et. al. identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease. In response to these data, the United States Food and Drug Administration (FDA) performed a patient level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9% - 23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2% - 14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16 - 2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices.

As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06 - 1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices.

The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

Additional information regarding clinical data:

In the ACOART I trial, the Kaplan Meier estimates freedom from all-cause mortality at 2, 3 and 5 years (within 5.5 years) are 91.9% (95%CI, 86.7% to 97.4%), 87.8% (95%CI, 81.6% to 94.5%), and 82.7% (95%CI, 75.6% to 90.5%), respectively, for the Acotec DCB treatment device and 93.9% (95%CI, 89.3% to 98.8%), 86.7% (95%CI, 80.2% to 93.7%) and 73.2% (95%CI, 64.0% to 83.6%), respectively, for the PTA control device.

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ACOART ORCHID DRUG ELUTING BALLOON CATHETER 0.035"							
Usable Catheter Length: 80 cm							
ORC03020AG	ORC35020AG	ORC04020AG	ORC45020AG	ORC05020AG	ORC55020AG	ORC06020AG	ORC07020AG
ORC03030AG	ORC35030AG	ORC04030AG	ORC45030AG	ORC05030AG	ORC55030AG	ORC06030AG	ORC07030AG
ORC03040AG	ORC35040AG	ORC04040AG	ORC45040AG	ORC05040AG	ORC55040AG	ORC06040AG	ORC07040AG
ORC03060AG	ORC35060AG	ORC04060AG	ORC45060AG	ORC05060AG	ORC55060AG	ORC06060AG	ORC07060AG
ORC03080AG	ORC35080AG	ORC04080AG	ORC45080AG	ORC05080AG	ORC55080AG	ORC06080AG	ORC07080AG
ORC03100AG	ORC35100AG	ORC04100AG	ORC45100AG	ORC05100AG	ORC55100AG	ORC06100AG	ORC07100AG
ORC03120AG	ORC35120AG	ORC04120AG	ORC45120AG	ORC05120AG	ORC55120AG	ORC06120AG	ORC07120AG
ORC03150AG	ORC35150AG	ORC04150AG	ORC45150AG	ORC05150AG	ORC55150AG	ORC06150AG	ORC07150AG
ORC03200AG	ORC35200AG	ORC04200AG	ORC45200AG	ORC05200AG	ORC55200AG	ORC06200AG	ORC07200AG
ORC03250AG	ORC35250AG	ORC04250AG	ORC45250AG	ORC05250AG	ORC55250AG	ORC06250AG	ORC07250AG
ORC03300AG	ORC35300AG	ORC04300AG	ORC45300AG	ORC05300AG	ORC55300AG	ORC06300AG	ORC07300AG

ORC08020AG	ORC09020AG	ORC10020AG	ORC12020AG
ORC08030AG	ORC09030AG	ORC10030AG	ORC12030AG
ORC08040AG	ORC09040AG	ORC10040AG	ORC12040AG
ORC08060AG	ORC09060AG	ORC10060AG	ORC12060AG
ORC08080AG	ORC09080AG	ORC10080AG	ORC12080AG
ORC08100AG	ORC09100AG	ORC10100AG	ORC12100AG
ORC08120AG	ORC09120AG	ORC10120AG	ORC12120AG

ACOART ORCHID DRUG ELUTING BALLOON CATHETER 0.035"							
Usable Catheter Length: 130 cm							
ORC03020BG	ORC35020BG	ORC04020BG	ORC45020BG	ORC05020BG	ORC55020BG	ORC06020BG	ORC07020BG
ORC03030BG	ORC35030BG	ORC04030BG	ORC45030BG	ORC05030BG	ORC55030BG	ORC06030BG	ORC07030BG
ORC03040BG	ORC35040BG	ORC04040BG	ORC45040BG	ORC05040BG	ORC55040BG	ORC06040BG	ORC07040BG
ORC03060BG	ORC35060BG	ORC04060BG	ORC45060BG	ORC05060BG	ORC55060BG	ORC06060BG	ORC07060BG
ORC03080BG	ORC35080BG	ORC04080BG	ORC45080BG	ORC05080BG	ORC55080BG	ORC06080BG	ORC07080BG
ORC03100BG	ORC35100BG	ORC04100BG	ORC45100BG	ORC05100BG	ORC55100BG	ORC06100BG	ORC07100BG
ORC03120BG	ORC35120BG	ORC04120BG	ORC45120BG	ORC05120BG	ORC55120BG	ORC06120BG	ORC07120BG
ORC03150BG	ORC35150BG	ORC04150BG	ORC45150BG	ORC05150BG	ORC55150BG	ORC06150BG	ORC07150BG
ORC03200BG	ORC35200BG	ORC04200BG	ORC45200BG	ORC05200BG	ORC55200BG	ORC06200BG	ORC07200BG
ORC03250BG	ORC35250BG	ORC04250BG	ORC45250BG	ORC05250BG	ORC55250BG	ORC06250BG	ORC07250BG
ORC03300BG	ORC35300BG	ORC04300BG	ORC45300BG	ORC05300BG	ORC55300BG	ORC06300BG	ORC07300BG

ORC08020BG	ORC09020BG	ORC10020BG	ORC12020BG
ORC08030BG	ORC09030BG	ORC10030BG	ORC12030BG
ORC08040BG	ORC09040BG	ORC10040BG	ORC12040BG
ORC08060BG	ORC09060BG	ORC10060BG	ORC12060BG
ORC08080BG	ORC09080BG	ORC10080BG	ORC12080BG
ORC08100BG	ORC09100BG	ORC10100BG	ORC12100BG
ORC08120BG	ORC09120BG	ORC10120BG	ORC12120BG

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
Acotec Scientific Co., Ltd
4th and 5th Floor, Building
No. 8, 16 Hongdabei Road,
BDA, 100176 Beijing, China
Tel : +86 10 67872107
Fax: +86 10 67866678

European Authorized
Representative:
MedDevConsult GmbH
Airport Center (Building C)
Flughafenstrasse 52a
22335 Hamburg
Germany

Publication Date : 2020 08 03